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PAGINĂ: 223 -

1. FOREWORD

The product complies with the essential requirements of 93/42/EEC Medical Device Directive and Transportable Pressure Equipment Directive 2010/35 EU.

The combination valve is designed according to EN ISO 10524-3 and EN ISO 10297 standards.

2. INTENDED USE

MediVital® Combination Valves are designed to be fitted to gas cylinders used for medical gases. These combination valves together with a gas cylinder form gas packages used either as gas supply point for medical devices (anaesthetic devices, ventilating devices, incubators etc.) or for direct gas supply to a patient's breathing mask or cannula.

GCE Combination Valves are intended to be used with the following medical gases:

- Oxygen
- Nitrous oxide
- Air for breathing
- Helium
- Carbon dioxide
- Mixtures of the gases listed above

MediVital® combination valve is not intended for use with flammable anaesthetics and substances.

2.1. PATIENT PROFILE

Medivital® supplies medical gas by means of the control flow head in flow range from 0 to 15 (25) l/min. Saturation of the medical gas must be prescribed by the doctor or an anaesthetic nurse who has got the medical education, assessed the health status of the patient and supplies him with required amount of the medical gas.

Health status - the assessment and decision about the use timing and dosage amount requirements of the patient is up to the doctor or medical personnel.

Medical personnel - nurse, rescuer, paramedic or carer who has obtained medical education.

2.2. USER PROFILE

Use in the hospital or ambulance:

- Education: personnel with medical education

Home care use:

Therapy is prescribed by a doctor or personnel who has got a medical education in the field of breathing support:

- Knowledge: training of personnel with medical education who has got an ability to read, train and apply knowledge from the instruction for use.

Training must be performed according to regulations valid in each country.

Important:

Patient must be trained by the doctor or person who has got an education in the life saving.

3. OPERATIONAL, TRANSPORT AND STORAGE SAFETY REQUIREMENTS

KEEP THE PRODUCT AND ITS ASSOCIATED EQUIPMENT AWAY FROM:

- Heat sources (fire, cigarettes,...)
- Flammable materials
- Oil or grease (take a great care in the use of hand creams)
- Water
- Dust

-  The product and its associated equipment must be prevented from tipping over, turning over or falling.
-  Always maintain oxygen cleanliness standards.
-  Only use the product and its associated equipment in well ventilated area.

Before the first use the product shall be kept in its original package. If removed from service (for transport, storage) GCE recommends using the original package (including inner packing materials).

National laws, rules and regulations for medical gases, accident prevention and environmental protection must be observed.

OPERATING CONDITIONS		STORAGE AND TRANSPORT CONDITIONS	
	-20* / +65 °C		-40 / +70 °C
	10 / 100 %		10 / 100 %
	600 / 1200 mbar		600 / 1200 mbar

*for inner tightness of the shut-off valve, during transport and storage of the combination valve mounted on a cylinder, the valid lower temperature limit is -40°C.

-  In case of combination valve storage at temperature below -20°C do not use the combination valve until its temperature reaches at least -20°C.
-  For the combination valves designed to be used with mixture of gases O₂+N₂O, the lowest operating temperature is +5°C. In normal use of the combination valve, frosting can appear on the combination valve surface, which is caused by the gas inside the combination valve when high pressure in the combination valve cooling when high pressure gas is being reduced to low pressure (Joule-Thomson effect). Check that all patient associated equipment connected to the combination valve is via a hose of at least 2 metres length.
-  O₂+N₂O mixtures are temperature sensitive. N₂O begins to separate out from the mixture if the temperature falls below about -6°C. A homogenous mixture is again obtained when the temperature has raised above 10°C and the cylinder was agitated. Before use, to ensure it is properly mixed, cylinders should be stored horizontally for 24 hours at a temperature above 10°C. If this is not practicable, before use the cylinders must be maintained at a temperature above 10°C for at least 2 hours and then completely inverted three times or placed in warm water at body temperature for 5 minutes and then completely inverted three times.

4. PERSONNEL INSTRUCTIONS

The Medical Devices Directive 93/42/EEC states that product provider must ensure that all personnel handling the product are provided with the operating instructions & performance data.

-  Do not use the product without properly familiarization of the product and its safe operation as defined in this Instruction for use. Ensure user is aware of particular information and knowledge required for the gas in use.



5. PRODUCT DESCRIPTION

Combination valve combines the function of shut-off valve of a high pressure gas cylinder and pressure regulator for the use with medical gases. Gas from the cylinder is first controlled by the main shut-off valve and then passed through the pressure regulator and finally delivered to the patient through the flow outlet or the pressure outlet. Outlet pressure is fixed by the manufacturer and each combination valve is provided with a low-pressure relief valve to protect against pressure regulator failure.

There are three basic alternatives:

- combination valve with quick-coupler outlet, the outlet pressure is constant, and flow is not controlled by the combination valve
- combination valve with outlet via calibrated nozzles, outlet flows is controlled by the combination valve and can be changed by the control head,
- the combination valve of both alternatives.

SHUT OFF VALVE STATUS INDICATOR:

Marking detail		
Inlet connection type	ISO V 30kg 17E-G5/8	Total weight allowance of the package
End of lifetime	2028/08 REF1 XXXXX	Filling port type
Reference number	REF 0718000 P _i =200 bar	Additional information
Serial number	SN 201308 XXXXX	Inlet pressure
Compliance with 2010/35 EU	CE XXXX TT0409	Compliance with 93/42/EEC

A – INLET STEM

The combination valve is fitted to gas cylinder by a threaded inlet stem. The inlet stem can be taper threaded or parallel threaded with different size depending on the cylinder size and material.

B – FILLING PORT

A filling port is provided for filling the gas cylinder at a filling station, it has no function for patient use. It includes a non-return valve (NRV). The NRV means that special filling adaptors are required to vent gas from the cylinder during the filling process (venting and/or vacuuming of cylinders).

C – INLET PRESSURE INDICATOR

Inlet pressure indicator is intended to indicate amount of gas in the gas cylinder. The pressure indicator is of an active type which means it indicates amount of gas in the gas cylinder whether the shut-off valve is opened or closed.

D – SHUT-OFF VALVE

The combination valve is provided with a shut-off valve to isolate the gas in the cylinder from the rest of combination valve functions. It must be opened during cylinder filling and patient's therapy. Part of the handwheel is an indicator showing open/close status of the shut-off valve.

The shut-off valve open/close status indicator is for guidance only. The shut-off valve may not be fully off when OFF status is showing. Fully closed status has to be checked by ensuring the shut-off valve is turned fully clockwise until a firm resistance to rotation is reached and gas supply to the outlets has stopped.

The shut-off valve must not be used in the "partly ON" status because even though gas is supplied to the outlets, the flow can be limited due to insufficiently opened shut-off valve.

E – RESIDUAL PRESSURE VALVE

Combination valve is equipped with a residual pressure valve with function to retain a minimum positive pressure in the gas cylinder to avoid contamination of the cylinder content by atmospheric air. During cylinder gas ventilation through the filling port the residual pressure valve is by passed.

F, G – FLOW CONTROL HEAD “F” AND FLOW OUTLET “G” (OPTION)

Combination valve can be delivered with a flow control head “F”. This function is used to supply adjustable gas flow rates (l/min) at atmospheric pressure directly to a patient through the flow outlet “G”, for instance through a cannula or a facemask. The flow outlet “G” is equipped with hose fitting (hose nipple) or a threaded type (for accessories to be connected via threaded connection).

Movement of the flow outlet ‘G’ is normal due to the method of fixing in the main body. It doesn't indicate a faulty flow outlet.

H – PRESSURE OUTLET OR QUICK COUPLER (OPTION)

The combination valve may be fitted with a pressure outlet. The pressure outlet is supplied with gas direct from the low-pressure part of the combination valve and it is fitted with a gas specific medical quick coupling connector also called “quick coupler”. The user can connect another piece of equipment to this outlet with a gas specific male probe (see appendix nr. 2). The quick coupler self seals when the male probe is disconnected. This outlet is for supplying gas at a controlled pressure to power medical devices, for instance medical ventilator.

I - PRESSURE RELIEVE VALVE OF LOW PRESSURE PART

Pressure relieve valve secures the low pressure part of combination valve and connected medical devices against over-pressure. If the gas pressure is decreased enough after the pressure relieve valve activation, it will closes itself.

J - EXCESS FLOW DEVICE OR DIP TUBE (OPTION)

Excess flow device ensures safe ventilation of gas from gas cylinder in case the combination valve is broken above inlet stem (e.g. cylinder fall). Dip tube does not have such function. Excess flow device and dip tube are to avoid contamination from cylinder entering the combination valve.

PRESSURE RELIEVE DEVICE OF HIGH PRESSURE PART - BURSTING DISC (OPTIONAL)

The high pressure relieve device is intended to protect the cylinder and high pressure part of combination valve against damage caused by increased cylinder pressure. If the pressure relieve device has been activated, it will not reseal and the combination valve must be taken out of service for repair (see Chapter 9).

Note: Colour of the product (especially guard, flow control head and shut-off valve) does not have to match the gas colour coding.

6. USE OF PRODUCT

6.1. OPERATIONS PERFORMED BY USER

6.1.1. BEFORE USE

VISUAL INSPECTION BEFORE USE:

- Check the combination valve for damage (incl. label and marking). If it shows signs of external damages, remove the product from service and suitably identify its status.
- Check the combination valve for contamination. If needed apply the cleaning procedure according to chapter 8.
- Check that the cylinder gas pressure indicator indicates sufficient pressure. If it indicates in the red zone, return the cylinder with combination valve back for filling.

LEAK-TIGHT AND FUNCTIONAL TEST BEFORE USE:

- Set flow control head (if any) on ZERO position - ensure the flow control head engages correctly.
- Slowly open shut-off valve handwheel (anticlockwise) until fully open - after approx. 1 turn (see indication in picture chapter 5).
- By listening check for leakage (leakage would be heard as characteristic hiss of flowing gas).
- Check that there is a gas flow at each flow control head set position in both clockwise and anticlockwise turning direction (for instance by sound or checking presence of bubbles in a humidifier).
- Close the shut-off valve (clockwise). Do not use excessive torque (max. closing torque is 5 Nm).
- Reset flow control head to ZERO position and ensure flow control head engages correctly.
- For combination valve fitted with pressure outlet, ensure it is in working condition by connecting and disconnecting quick-coupler probe.

6.1.2. USE OF COMBINATION VALVE

6.1.2.1. Use of combination valve flow outlet and setting of flow:

- Ensure the flow control head is in position "0".
- Connect accessories to the flow outlet.
- Slowly open shut-off valve handwheel (anticlockwise) until fully open - after approx. 1 turn (see indication in picture chapter 5).
- Set flow head to required flow position. Ensure flow control head engages correctly.

Before connecting any accessory to the flow outlet make sure that the patient is not connected.

Sudden opening could result in a danger of fire or explosion arising from oxygen pressure shocks. Insufficient opening of the shut-off valve could reduce actual flow delivered.

Always ensure that the flow control head has been correctly set and not placed between two settings. Flow control head placed between settings will not deliver correct flow of medical gas.

 Common variants of flow control head can have an “end stop” in between the maximum flow position and the zero position. Do not try to apply excessive force on the flow control head when it stops in the maximum flow position (during clockwise rotation) or in zero position (anti-clockwise rotation).

 Medical gas flow rate must be prescribed by a doctor.

6.1.2.2. Use of combination valve pressure outlet

- Ensure the flow control head is in position “0” (if included).
- Slowly open shut-off valve handwheel (anticlockwise) until fully open - after approx. 1 turn (see indication in picture chapter 5).
- Connect accessory to the pressure outlet if not already connected.

 Before connecting any accessory to the pressure outlet make sure that the patient is not connected and the accessory outlet is secured.

 If pressure outlet is to be connected to a medical device that requires high gas flow (for instance pulmonary ventilator that requires gas flow 100 l/min at the minimum pressure 2.8 bar), compare the required flow of connecting medical device with pressure and flow characteristics of the combination valve stated in appendix Nr. 1. To assure sufficient performance (pressure and flow characteristics of the combination valve) the medical device should not be used if pressure indicator enters the red zone.

 Sudden opening could result in a danger of fire or explosion arising from oxygen pressure shocks. Insufficient opening of the shut-off valve could reduce actual flow delivered.

 If a pressure outlet as well as a flow outlet are part of the combination valve do not use them simultaneously, especially if pressure in the cylinder is below 50 bar, in could adversely affect the outlet parameters of the combination valve.

6.1.3. AFTER USE

- Close shut-off valve (clockwise). Do not use excessive torque(max. closing torque is 5 Nm).
- Vent pressure from the connected devices.
- Disconnect all connected devices from user outlets.
- Set flow control head on "0" (if included).

7. ACCESSORIES

ACCESSORIES CONNECTABLE TO FLOW OUTLET:

- hose connected with mask, cannula or humidifier.

ACCESSORIES CONNECTABLE TO PRESSURE OUTLET:

- low pressure hose (working pressure >10 bar), flowmeters, Venturi suction ejectors, lung ventilators.

OTHER USER ACCESSORIES:

- bed hanger, humidifier holder.

ACCESSORIES FOR FILLING STATIONS:

- filling adaptor.

 Before connecting any accessory or medical device to the combination valve, always check that they are fully compatible with connection features & performances of the product.

8. CLEANING

Remove dirt with a soft cloth dampened in oil free soap water & rinsed with clean water. Disinfection can be carried out with an alcohol-based solution (with dampened wipes). If other cleaning solutions are used, check that they are not abrasive and they are compatible with the product materials (including labels) and gas (convenient cleaning solution - i.e. Meliseptol).

-  Do not use cleaning solutions containing ammonia!
-  Do not expose to water or any other liquid.
-  Do not expose to high temperature (such as autoclave).
-  To apply the cleaning solution do not spray it as the spray may enter into the inner parts of combination valve and cause contamination or damage.
-  Do not use pressure wash as it could damage or contaminate the combination valve.
-  If the inner parts of the combination valve have been contaminated do not continue to use the combination valve under any circumstances. It must be withdrawn from service.

9. SERVICE, PRODUCT LIFE TIME AND MAINTENANCE

9.1. PRODUCT LIFE TIME

SERIAL NUMBER AND MANUFACTURE DATE

An eleven-figure serial number stamped on combination valve body consists of the data as follows:

YYYYMM XXXXX

YYYY: year of manufacture

MM: month of manufacture

XXXXX: sequence number of product

For example: Serial number 201303 00521 indicates the combination valve manufactured in 2013, in March, with sequence number 521.

PRODUCT LIFE TIME AND WASTE MANAGEMENT

End of life time is marked on the combination valve body this way:

YYYYMM

YYYY: year of lifetime end

MM: month of lifetime end

Maximum life time of this product is 15 years from date of manufacturing.

At the end of the product's life time, the product must be withdrawn from service. The provider of the device shall prevent the reuse of the product and handle the product in compliance with "Directive of European Parliament and Council 2008/98/EC on waste".

In accordance to Article 33 of REACH GCE, s.r.o. as responsible manufacturer shall inform all customers if materials containing 0.1% or more of substances included in the list of Substance of Very High Concern (SVHC). The most commonly used brass alloys used for valve bodies and other brass components contain 2-3% of lead (Pb), EC no. 231-468-6, CAS no. 7439-92-1. The lead will not be released to the gas or surrounding environment during normal use. After end of life the product shall be scrapped by an authorized metal recycler to ensure efficient material handling with minimal impact to environment and health. To date we have no information that indicates that other materials containing SVHC of concentrations exceeding 0.1% are included in any GCE product.

9.2. REPAIR AND SERVICE OPERATIONS

Repairs and service can be only done by a GCE certified person who also holds all necessary certificates in compliance with national standards for mounting and repair of dedicated gas devices. For information about a service in your area please contact GCE or distributor of GCE product. Usually the combination valves can be repaired while fitted on the cylinder.

Repair that do not need to be done by certified personnel include exchange of the below mentioned components:

- guard,
- labels,
- protective covers and separable hose adaptors.

 **All labels on the equipment must be kept in good, legible condition by the provider and the user during the entire product life time.**

 **Use only genuine GCE components.**

Any product sent back to GCE (or GCE authorised centre) for repair or maintenance shall be properly packaged to prevent contamination or damage during storage, transportation and handling.

For product to be repaired, the fault short description or some reference to a claim nr. should be indicated.

10. GLOSSARY

	Consult instruction for use		Suitable for Hospital care use
	Caution!		Suitable for Home care use
	Keep away from heat and flammable materials		Suitable for Emergency care use
	Keep away from oil and grease!		Fragile, handle with care
	Keep dry!		Weight of product
	Humidity limit		Atmospheric pressure limitation
	Temperature limit		Use by date
	Inlet parameter		Serial number
	Outlet parameter		Catalogue number
P_1	Inlet pressure		Batch code
P_2	Outlet pressure		Customer number

	Manufacturer		Take the equipment back for recycling. Do not dispose equipment into unsorted municipal waste.
	Date of manufacture		

EN

11. WARRANTY

The Standard Warranty period is two years from date of receipt by the GCE Customer (or if this is not known 2 years from time of the product manufacture shown on the product).

The standard warranty is only valid for products handled according to Instruction for use (IFU) and general industry good practice and standards.

APPENDIX:

Nr 1 – Technical and performance data

Nr 2 – Quick-coupler features and connecting/disconnecting procedure

Nr 3 – Post Filling Checks

Nr 4 – Valve assembly and filling instruction

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© GCE, s.r.o.